

Morinaga Food Allergen Lateral Flow  
Rapid Test Pro II for Egg

Rev. Jan. 2022

Lot number from #PE001 -

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result
<b>Cereals</b>	Rice powder	-	<b>Sea foods</b>	Salmon	-	<b>Coffee·Cocoa·Tea</b>	Coffee bean	-
	Polished rice	-		Mackerel	-		Cocoa powder	-
	Sprouted brown rice	-		Squid	-		Green tea	-
	Glutinous rice	-		Octopus	-	Roasted laver	-	
	Buckwheat	-		Black tiger prawn	-	Brown seaweed	-	
	Japanese barnyard millet	-		Snow crab	-	Kelp 200-fold extract	+	
	Foxtail millet	-		Scallop adductor	-	<b>Spices</b>	Cumin	-
	Common millet	-		Pollack roe	-		Poppy seeds	-
	Wheat	-		Salmon roe	-		Black pepper	-
	Rye	-		Beef	-		White pepper	-
	Barley	-	Chicken	-	Red pepper		-	
	Malt	-	Pork	-	Wasabi (green horseradish)		-	
	Oat	-	Chicken egg	-	Ginger		-	
	Adlay※	-	Processed chicken egg	-	Curry powder (spice mix)		-	
	Corn	-	Dried chicken egg	-	Cinnamon		+	
	Corn flour	-	<b>Dairy products</b>	Milk	-		Coriander	-
	Amaranth	-		Skim milk	-		Turmeric	-
	Quinoa	-		Tomato	-		Polysaccharide Thickener	Low methoxyl pectin 50-fold extract
<b>Tubers·Starch</b>	Sago palm powder	-	<b>Green stuff</b>	Onion	-	<b>Others</b>	Guar gum 1000-fold extract	-
	Yam	-		Spinach	-		Gelatin	-
	Potato starch	-	<b>Fungus</b>	Mushroom	-			
<b>Beans</b>	Soybean	-	<b>Fruits</b>	Orange	-			
	Red bean (azuki)	-		Kiwi	-			
	Common bean	-		Apple	-			
	Concord pole	-		Peach	-			
	Broad bean	-		Banana	-			
Peanut	-	Coconut milk		-				
<b>Seeds</b>	Roasted almond	-						
	Roasted cashew nut	-						
	Roasted macadamia nut	-						
	Roasted pistachio nut	-						
	Roasted walnut	-						
	White sesame	-						
Black sesame	-							

Test Result	
-	Negative
+	Positive (reaction positive)

N.T. : Not Tested

※The Test Solution may be difficult to flow on a test stick.

These samples were tested by following the instruction of Rapid Test Pro II for Egg (using Recommended Extraction Method), except Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.