

Morinaga Food Allergen Lateral Flow
Rapid Test Pro II for Egg

Rev. Sep. 2022

Lot number from #PE023 -

Group	Ingredient	Test Result	Group	Ingredient	Test Result	Group	Ingredient	Test Result	
Cereals	Rice powder	-	Seeds	Peanut	-	Fruits	Orange	-	
	Polished rice	-		Roasted almond	-		Kiwi	-	
	Sprouted brown rice	-		Roasted cashew nut	-		Apple	-	
	Glutinous rice	-		Roasted macadamia nut	-		Peach	-	
	Buckwheat	-		Roasted pistachio nut	-		Banana	-	
	Japanese barnyard millet	-		Roasted walnut	-		Coconut milk	-	
	Foxtail millet	-		Roasted pecan nut	-	Coffee bean	-		
	Common millet	-		White sesame	-	Cocoa powder	-		
	Wheat	-		Black sesame	-	Green tea	-		
	Rye	-		Sea foods	Salmon	-	Roasted laver	-	
	Barley	-	Mackerel	-	Hijiki (boiled)	-			
	Malt	-	Squid	-	Brown seaweed (raw)	- *			
	Oat	-	Octopus	-	Kelp 200-fold extract	-			
	Adlay	-	Black tiger prawn	-	Cumin	-			
	Corn	-	Snow crab	-	Poppy seeds	-			
	Corn flour	-	Scallop adductor	-	Black pepper	-			
Amaranth	-	Pollack roe	-	White pepper	-				
Quinoa	-	Salmon roe	-	Red pepper	-				
Tubers· Starch	Sago palm powder	-	Meats	Beef	-	Spices	Ginger	-	
	Yam	-		Chicken	-		Cinnamon	-	
	Potato starch	-		Pork	-		Coriander	-	
Beans	Soybean	-	Eggs	Chicken egg	+		Polysaccharide Thickener	Turmeric	-
	Red bean (azuki)	-		Processed chicken egg	+			Low methoxyl pectin 50-fold extract	-
	Common bean	+		Dried chicken egg	+		Guar gum 1000-fold extract	-	
	Concord pole	-	Dairy products	Milk	-		Others	Gelatin	-
	Broad bean	-		Skim milk	-				
		Green stuff	Tomato	-					
			Onion	-					
			Spinach	-					
		Fungus	Mushroom	-					

Test Result	
-	Negative
+	Positive (reaction positive)

*The Test Solution may be difficult to flow on a test stick.

These samples were tested by following the instruction of Rapid Test Pro II (using Recommended Extraction Method), except Kelp 200-fold extract, Polysaccharide and Thickener samples.

These results can be changed by the difference of kit's lot number or the sample's species or production area.